

Just 3 Minutes to More Cost-Effective Lead Testing


*Magellan Diagnostics offers the
fastest, easiest way to test for lead.*

With the widest range of lead
testing systems available, we have
a solution that's right for you.



Blood Lead Testing In-House

Which system is right for you?



	LeadCare II	LeadCare Plus*	LeadCare Ultra
Typical setting	Physician office	Small / mid-size lab	Large lab
CLIA complexity	Waived	Moderate	Moderate
Sample throughput	n/a	15–20 / hour	Up to 90 / hour
Analytical range	3.3–65 µg/dL	1.9–65 µg/dL	1.9–65 µg/dL
Confirmation by reference lab required?	Yes	No	No
Sample stability	n/a	Up to 72 hours post draw	Up to 72 hours post draw
Sample type	Capillary or venous whole blood (EDTA or heparin)	Capillary or venous whole blood (EDTA or heparin)	Capillary or venous whole blood (EDTA or heparin)
Sample Collection	Capillary tubes (included in kit)	Microtainer or Vacutainer (not included)	Microtainer or Vacutainer (not included)
Number of channels	1	1	6
Tests per kit	48	96	192
Connectivity	No	Optional data management system (HL7 2.5.1)	Yes

* Pending FDA 510(k) clearance.

The Most Trusted Name In Lead Testing

The LeadCare platform has been used in clinical labs for nearly twenty years. We are proud to provide a family of FDA-cleared systems offering the fastest, simplest way to perform clinical blood lead testing.